

CLINICAL RESULTS

For over two decades, InMode has revolutionized the medical aesthetic industry with state-of-the-art technologies. InMode now has over 50 peer-reviewed and textbook publications, leading the industry in clinical excellence.



Major Beauty Awards 2020



WWW.INMODEMD.COM

the new



MORPHEUS PRO

by INMODE



DISCOVER NEW DEPTHS

Morpheus8 is InMode's newest epidermal resurfacing and sub-dermal adipose remodeling device (SARD). Morpheus8 is the first and only programmable, full-body, fractional technology to morph the aging face and contour the body.

The triple action of sub-necrotic bulk heating, coagulation of fat and connective tissue, and ablation results in improvement in skin surface irregularities, skin tone, laxity, as well as treating fine lines, wrinkles and indented scars.

- Smooths and tightens fine lines and wrinkles
- Treats depressed scars (acne)
- Improves active acne
- Smooths skin irregularities
- Improves skin tone
- Lifts and contours face, neck and body areas



FULL BODY FRACTIONAL
REMODELING

MORPHEUS8

BENEFITS

- Performs deep and safe fractional treatment with programmable penetration from 1-7mm, depending on the configuration
- Polymer coated micropins assure virtually no thermal damage to epidermis
- Safely treats skin type VI with little risk of PIH
- Disposable tip prevents cross-contamination
- Advanced lattice design of return electrodes creates a uniform effect

The Morpheus8 is a fractional device for treatment of facial surface irregularities, such as fine lines & wrinkles and acne scars, together with body remodeling for treatment of cellulite and improving general skin tone. The range of Morpheus8 tip configurations, with programmable penetration depth and energy delivery for targeted treatments, is ideal for a complete facial and body remodeling regimen.

When used together with InMode's BodyTite, FaceTite and AccuTite handpieces for radio frequency-assisted lipolysis (RFAL), the Morpheus8 triple action of fat coagulation, connective tissue contraction and sub-necrotic heating in the entire subdermis and dermis results in the perfect solution for patients desiring a complete facial and body reshaping regimen, with immediate results, and continued improvement over time.

PRIME

12 gold-coated micropins



MORPHEUS8

24 gold-coated micropins



RESURFACING

24 gold-coated micropins



BODY

40 gold-coated micropins



	Prime	Morpheus8	Morpheus8 Resurfacing	Body
Treatment zones	Subdermal and dermal remodeling of focal areas through fractional coagulation and sub-necrotic bulk heating	Subdermal and dermal remodeling for face, neck and body areas	Epidermal resurfacing of the superficial skin layer	Dermal and subdermal remodeling of large body areas such as abdomen, thighs and buttocks through fractional coagulation and bulk heating.
No. pins	12 coated micropins	24 coated micropins	24 coated micropins	40 coated micropins
Pin length	1 - 4 mm, programmable	1 - 4 mm programmable	0.5mm fixed	2 - 7 mm, programmable
Depth of treatment	1 - 4 mm programmable	Facial & bony areas - 2mm Soft facial areas - 3mm Body - 4mm	0.5mm	2 - 7 mm, programmable
Coagulation depth	Up to 4mm	Up to 4 mm	0.5 mm	Up to 7 mm
Heating depth	Up to 1mm additional heat profile	Up to 1mm additional heat profile	Up to 0.75mm additional heat profile	Up to 1mm additional heat profile


Introducing MORPHEUS8 BODY

Morpheus8 Body provides deep 8mm subdermal adipose tissue remodeling for larger body areas such as thighs, buttocks and abdomen at the deepest levels of RF energy available. The Morpheus8 Body tip consists of a matrix of 40 gold-coated microneedles, with a maximum depth penetration of 7mm and an additional heat profile of 1 mm.



Three step RF Vulcanization Process:


- CHARGE - Powerful RF energy burst delivered deep into the FSN causing re-orientation of collagen fibers
- MOLD - RF energy dissipation leads to block heating resulting in collagen remodeling and buttressing of FSN
- WEAVE - RF energy delivered to the dermal level causing pulling and synching of fibers into a weave



Morpheus8 Body Burst Customizable multi-pulse settings

Specifications	Cycle Setting	Mult-Level Penetration Cycle
	7mm	7mm, 5mm, 3mm
	6mm	6mm, 4mm, 2mm
	5mm	5mm, 3mm
	4mm	4mm, 2mm

DEPTH OF PENETRATION



FORMA



The first auto-adjusting thermal skin treatment for collagen remodeling and skin contraction for all areas of the body

Forma is non-invasive, sub-dermal heating applicator for collagen remodeling. RF current flows uniformly between the electrodes to provide a comfortable thermal experience. Built-in temperature controls allow for the attainment of clinically proven optimal temperatures, and the exclusive capacity for prolonged exposure at therapeutic temperatures.

Specifications	RF Output Power	Up to 65 W
	Output Frequency	1 MHz
	Temperature Cut-Off	Real time temperature measurements set by operator up to 43°C